

14-3-3 sigma Antibody (Biotin Conjugate) (R33530BTN)

Catalog No.	Formulation	Size
R33530BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Pig, Rat
Format	Biotin Conjugate
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2810
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This 14-3-3 sigma antibody is available for research use only.



Western blot testing of mouse skin lysate with biotinylated 14-3-3 sigma antibody at 0.1ug/ml. Predicted molecular weight: ~28kDa.

Description

Additional name(s) for this target protein: Stratifin; SFN

Application Notes

Optimal dilution of the 14-3-3 sigma antibody should be determined by the researcher.

Immunogen

Amino acids DLHTLSEDSYKDST were used as the immunogen for this 14-3-3 sigma antibody.

Storage

Aliquot and store the 14-3-3 sigma antibody at -20oC.