


14-3-3 zeta/delta Antibody / YWHAZ [clone 6H7] (RQ5643)

Catalog No.	Formulation	Size
RQ5643	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat, Monkey
Format	Antigen affinity purified
Clonality	Monoclonal
Isotype	Mouse IgG1
Clone Name	6H7
Purity	Affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose and 0.025% sodium azide
UniProt	P63104
Applications	Western Blot : 0.5-1ug/ml
Limitations	This 14-3-3 zeta/delta antibody is available for research use only.



Western blot testing of 1) human HeLa, 2) human A549, 3) monkey COS-7, 4) human Raji, 5) human Caco-2, 6) human Jurkat, 7) mouse brain and 8) rat brain lysate with 14-3-3 zeta/delta antibody. Predicted molecular weight ~28 kDa.

Description

14-3-3 protein zeta/delta is a protein that in humans is encoded by the YWHAZ gene on chromosome 8. This gene product belongs to the 14-3-3 family of proteins which mediate signal transduction by binding to phosphoserine-containing proteins. This highly conserved protein family is found in both plants and mammals, and this protein is 99% identical to the mouse, rat and sheep orthologs. The encoded protein interacts with IRS1 protein, suggesting a role in regulating insulin sensitivity. Several transcript variants that differ in the 5' UTR but that encode the same protein have been identified for this gene.

Application Notes

Optimal dilution of the 14-3-3 zeta/delta antibody should be determined by the researcher.

Immunogen

Amino acids LLEKFLIPNASQAESKVFYLMKGDYYRYLAEVAAGDDKKGIVDQ from the human protein were used as the immunogen for the 14-3-3 zeta/delta antibody.

Storage

After reconstitution, the 14-3-3 zeta/delta antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.