


5HT3B receptor Antibody (R35996)

Catalog No.	Formulation	Size
R35996-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9177
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This 5HT3B receptor antibody is available for research use only.



Western blot testing of human cerebral cortex lysate with 5HT3B receptor antibody at 1ug/ml. Predicted molecular weight (human): ~50kDa.

Description

Additional name(s) for this target protein: 5-hydroxytryptamine (serotonin) receptor 3B; HTR3B

Application Notes

Optimal dilution of the Serotonin Receptor 3B / HTR3B antibody should be determined by the researcher.

Immunogen

Amino acids SNYLQTQDQTDQQE were used as the immunogen for this 5HT3B receptor antibody.

Storage

Aliquot and store the 5HT3B receptor antibody at -20oC.