

A1BG Antibody (F50307)


Catalog No.	Formulation	Size
F50307-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F50307-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Purified
UniProt	P04217
Applications	Western Blot : 1:1000 IHC (Paraffin) : 1:50-1:100
Limitations	This A1BG antibody is available for research use only.



Western blot analysis of A1BG antibody and human T47D lysate. Predicted molecular weight ~55 kDa.



A1BG antibody IHC analysis in formalin fixed and paraffin embedded human breast carcinoma tissue.

Description

Alpha-1B-glycoprotein is a plasma glycoprotein of unknown function. The protein shows sequence similarity to the variable regions of some immunoglobulin supergene family member proteins.

Application Notes

Titration of the A1BG antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 421-449 from the human protein was used as the immunogen for this A1BG antibody.

Storage

Aliquot the A1BG antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.