


## ABCC4 Antibody / MRP4 (R33289)

Catalog No.	Formulation	Size
R33289-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10257
Applications	Immunohistochemistry (FFPE) : 2-4ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This ABCC4 antibody is available for research use only.



IHC testing of FFPE human prostate with ABCC4 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

### Description

Additional name(s) for this target protein: ATP-binding cassette, sub-family C

### Application Notes

Optimal dilution of the ABCC4 antibody should be determined by the researcher.

1. This ABCC4 antibody is specific for isoform 1 only.

## Immunogen

Amino acids CNGQPSTLTIFETAL were used as the immunogen for this ABCC4 antibody.

## Storage

Aliquot and store the ABCC4 antibody at -20oC.