


ABCD2 Antibody (F42334)


Catalog No.	Formulation	Size
F42334-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F42334-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

Bulk quote request


Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q9UBJ2
Localization	Cytoplasmic
Applications	Western Blot : 1:1000 IHC (Paraffin) : 1:10-1:50 Flow Cytometry : 1:10-1:50
Limitations	This ABCD2 antibody is available for research use only.



ABCD2 antibody immunohistochemistry analysis in formalin fixed and paraffin embedded human heart tissue.



ABCD2 antibody flow cytometric analysis of K562 cells (right histogram) compared to a [negative control](http://www.ncbi.nlm.nih.gov/pmc/articles/PMC1001001/) (left histogram). FITC-conjugated donkey-anti-rabbit secondary Ab was used for the analysis.



ABCD2 antibody western blot analysis in K562 lysate. Predicted molecular weight: ~83 kDa.

Description

The protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of the ALD subfamily, which is involved in peroxisomal import of fatty acids and/or fatty acyl-CoAs in the organelle. All known peroxisomal ABC transporters are half transporters which require a partner half transporter molecule to form a functional homodimeric or heterodimeric transporter. The function of this peroxisomal membrane protein is unknown; however this protein is speculated to function as a dimerization partner of ABCD1 and/or other peroxisomal ABC transporters. Mutations in this gene have been observed in patients with adrenoleukodystrophy, a severe demyelinating disease. This gene has been identified as a candidate for a modifier gene, accounting for the extreme variation among adrenoleukodystrophy phenotypes. This gene is also a candidate for a complement group of Zellweger syndrome, a genetically heterogeneous disorder of peroxisomal biogenesis.

Application Notes

Titration of the ABCD2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 553-582 from the human protein was used as the immunogen for this ABCD2 antibody.

Storage

Aliquot the ABCD2 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.