


ABTB2 Antibody (F42258)

Catalog No.	Formulation	Size
F42258-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F42258-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q8N961
Applications	Western Blot : 1:1000
Limitations	This ABTB2 antibody is available for research use only.



Description

Ankyrin repeat and BTB/POZ domain-containing protein 2 may be involved in the initiation of hepatocyte growth (By similarity).

Application Notes

Titration of the ABTB2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 681-710 from the human protein was used as the immunogen for this ABTB2 antibody.

Storage

Aliquot the ABTB2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.