


ACADM Antibody / Mcad (R35671)


Catalog No.	Formulation	Size
R35671-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	34
Applications	Western Blot : 0.01-0.1ug/ml Immunohistochemistry (FFPE) : 3-5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This ACADM antibody is available for research use only.



Western blot testing of human heart lysate with ACADM antibody at 0.05ug/ml. Predicted molecular weight: ~46kDa.



IHC staining of FFPE human liver with ACADM antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Acyl-Coenzyme A dehydrogenase medium chain

Application Notes

Optimal dilution of the ACADM antibody should be determined by the researcher.

Immunogen

Amino acids RLIVAREHIDKYKN were used as the immunogen for this ACADM antibody.

Storage

Aliquot and store the ACADM antibody at -20oC.