


## ACK1 Antibody / Activated CDC42 kinase 1 (F50698)


Catalog No.	Formulation	Size
F50698-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F50698-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human, Mouse
<b>Predicted Reactivity</b>	Bovine
<b>Format</b>	Purified
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Purified
<b>UniProt</b>	Q07912
<b>Localization</b>	Cytoplasmic
<b>Applications</b>	Western Blot : 1:500-1000 IHC (Paraffin) : 1:50-1:100
<b>Limitations</b>	This ACK1 antibody is available for research use only.



Western blot testing of human brain lysate with ACK1 antibody at 1:1000. Predicted molecular weight ~114 kDa.



IHC analysis of FFPE human hepatocarcinoma tissue stained with the ACK1 antibody

## Description

Activated CDC42 kinase 1 is a tyrosine kinase that binds Cdc42Hs in its GTP-bound form and inhibits both the intrinsic and GTPase-activating protein (GAP)-stimulated GTPase activity of Cdc42Hs. This binding is mediated by a unique sequence of 47 amino acids C-terminal to an SH3 domain. The protein may be involved in a regulatory mechanism that sustains the GTP-bound active form of Cdc42Hs and which is directly linked to a tyrosine phosphorylation signal transduction pathway.

## Application Notes

Titration of the ACK1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## Immunogen

A portion of amino acids 934-964 from the human protein was used as the immunogen for this ACK1 antibody.

## Storage

Aliquot the ACK1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.