

Aconitase 1 Antibody (R35495)

Catalog No.	Formulation	Size
R35495-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Human, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	48
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This Aconitase 1 antibody is available for research use only.



Western blot testing of human liver lysate with Aconitase 1 antibody at 0.1ug/ml.
Predicted molecular weight: ~98kDa.

Description

Additional name(s) for this target protein: Aconitate hydratase 1; ACO1

Application Notes

Optimal dilution of the Aconitase 1 antibody should be determined by the researcher.

Immunogen

Amino acids QVDFNRRADSLQKNQ were used as the immunogen for this Aconitase 1 antibody.

Storage

Aliquot and store the Aconitase 1 antibody at -20oC.