


ACOX2 Antibody (R35514)

Catalog No.	Formulation	Size
R35514-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Rabbit
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8309
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This ACOX2 antibody is available for research use only.



Western blot testing of human liver lysate with ACOX2 antibody at 0.5ug/ml with [B] and without [A] blocking/immunizing peptide. Predicted molecular weight: ~75kDa.

Description

Additional name(s) for this target protein: Acyl-Coenzyme A oxidase 2

Application Notes

Optimal dilution of the ACOX2 antibody should be determined by the researcher.

Immunogen

Amino acids QWAQKSPTNTQENP were used as the immunogen for this ACOX2 antibody.

Storage

Aliquot and store the ACOX2 antibody at -20oC.