


ACPP Antibody / Prostatic Acid Phosphatase / PAP (F50540)


Catalog No.	Formulation	Size
F50540-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F50540-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Purified
UniProt	P15309
Localization	Cytoplasmic
Applications	Western Blot : 1:1:500-1000 IHC (Paraffin) : 1:10-1:50
Limitations	This ACPP antibody is available for research use only.



IHC analysis of FFPE human lung carcinoma tissue stained with ACPP antibody



Western blot analysis of ACPP antibody and mouse lung tissue lysate. Expected molecular weight: 45-50 kDa.

Description

Prostatic acid phosphatase is an enzyme which catalyzes the conversion of orthophosphoric monoester to alcohol and orthophosphate. It is synthesized under androgen regulation and is secreted by the epithelial cells of the prostate gland.

Application Notes

Titration of the ACPD antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 313-341 from the human protein was used as the immunogen for this ACPD antibody.

Storage

Aliquot the ACPD antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.