


ADAM19 Antibody (R34477)


Catalog No.	Formulation	Size
R34477-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Pig, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8728
Applications	Western Blot : 0.1-0.3ug/ml IHC (FPPE) : 5ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This ADAM19 antibody is available for research use only.



Western blot testing of spleen lysate with ADAM19 antibody at 0.1ug/ml. Predicted molecular weight: 105-125 kDa depending on level of glycosylation.



IHC testing of FFPE human spleen tissue with ADAM19 antibody at 5ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

Description

Additional name(s) for this target protein: ADAM metallopeptidase domain 19, MADDAM

Application Notes

Optimal dilution of the ADAM19 antibody should be determined by the researcher.

Immunogen

Amino acids KIQCQSSEARPLESN were used as the immunogen for this ADAM19 antibody.

Storage

Aliquot and store the ADAM19 antibody at -20oC.