

Adenosine A2a Receptor Antibody (R34176)

Catalog No.	Formulation	Size
R34176-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Rhesus Monkey
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	135
Applications	Western Blot : 1-3ug/ml (1) ELISA (peptide) LOD : 1:8000
Limitations	This Adenosine A2a Receptor antibody is available for research use only.



Description

Additional name(s) for this target protein: RDC8, hA2aR; ADORA2A

Application Notes

Optimal dilution of the Adenosine A2a Receptor antibody should be determined by the researcher.

1. Reported to work in primary Rhesus microglia lysate by western blot.

Immunogen

Amino acids HELKGVCPPEPGLDD were used as the immunogen for this Adenosine A2a Receptor antibody.

Storage

Aliquot and store the Adenosine A2a Receptor antibody at -20oC.