


Adipose triglyceride lipase Antibody (R34960)


Catalog No.	Formulation	Size
R34960-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	57104
Applications	Immunofluorescence : 10ug/ml Flow Cytometry : 1ug/million cells in 0.1ml ELISA (peptide) LOD : 1:8000
Limitations	This Adipose triglyceride lipase antibody is available for research use only.



IF/ICC staining of fixed and permeabilized human HEK293 cells with Adipose triglyceride lipase antibody (green) at 10ug/ml and DAPI nuclear stain (blue).



FACS testing of fixed and permeabilized human HEK293 cells with Adipose triglyceride lipase antibody (blue) at 1ug/million cells/0.1ml and naive goat Ig (black).

Description

Additional name(s) for this target protein: ATGL; PNPLA2

Application Notes

Optimal dilution of the Adipose triglyceride lipase antibody should be determined by the researcher.

Immunogen

Amino acids PARPHGPEDKDQ were used as the immunogen for this Adipose triglyceride lipase antibody.

Storage

Aliquot and store the Adipose triglyceride lipase antibody at -20oC.