


AES Antibody (R35917)

Catalog No.	Formulation	Size
R35917-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	166
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This AES antibody is available for research use only.



Western blot testing of human cerebral cortex lysate with AES antibody at 0.3ug/ml.
Predicted molecular weight: ~22 kDa (AES-1), ~29 kDa (AES-2).

Description

Additional name(s) for this target protein: Amino-terminal enhancer of split, AES-1, AES-2

Application Notes

Optimal dilution of the AES antibody should be determined by the researcher.

Immunogen

Amino acids SHLPQQLKFTTSD were used as the immunogen for this AES antibody.

Storage

Aliquot and store the AES antibody at -20oC.