


AGRP Antibody (R33323)

Catalog No.	Formulation	Size
R33323-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	181
Applications	Immunohistochemistry (Frozen) : 0.05-0.1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This AGRP antibody is available for research use only.



IHC testing of PFA-perfused cryosection of human hypothalamus and AGRP antibody at 0.05ug/ml. HIER: pH6 citrate buffer at 80C for 30min, HRP-staining with Ni-DAB after biotin-SP-anti-goat amplification.

Description

Additional name(s) for this target protein: Agouti related protein

Application Notes

Optimal dilution of the AGRP antibody should be determined by the researcher.

Immunogen

Amino acids CRKLG TAMNPCSRT were used as the immunogen for this AGRP antibody.

Storage

Aliquot and store the AGRP antibody at -20oC.