


AHSG Antibody (R35387)

Catalog No.	Formulation	Size
R35387-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	197
Applications	Western Blot : 0.2-0.5ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This AHSG antibody is available for research use only.



Western blot testing of human albumin-depleted serum with AHSG antibody at 0.5ug/ml.
Predicted molecular weight: ~39/45-59kDa (unmodified/glycosylated).

Description

Additional name(s) for this target protein: Alpha-2-HS-glycoprotein, fetuin-A

Application Notes

Optimal dilution of the AHSG antibody should be determined by the researcher.

Immunogen

Amino acids QPEGANEAVPTP were used as the immunogen for this AHSG antibody.

Storage

Aliquot and store the AHSG antibody at -20oC.