


AID Antibody (R35749)

Catalog No.	Formulation	Size
R35749-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	57379
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This AID antibody is available for research use only.



Western blot testing of human tonsil lysate with AID antibody at 0.1ug/ml. Predicted molecular weight ~24kDa.

Description

Additional name(s) for this target protein: Activation-induced deaminase, AICDA, activation-induced cytidine deaminase; AICDA

Application Notes

Optimal dilution of the AID antibody should be determined by the researcher.

Immunogen

Amino acids RRKFLYQFKNVR were used as the immunogen for this AID antibody.

Storage

Aliquot and store the AID antibody at -20oC.