


AIF-1 Antibody / Allograft inflammatory factor 1 / IBA1 [clone AFGH-1] (RQ8894)

Catalog No.	Formulation	Size
RQ8894	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	AFGH-1
Purity	Affinity chromatography
UniProt	P55008
Localization	Cytoplasm
Applications	Western Blot : 1:500 Immunohistochemistry (FFPE) : 1:50
Limitations	This AIF-1 antibody is available for research use only.



IHC staining of FFPE human tonsil tissue with AIF-1 antibody, HRP-secondary and DAB substrate. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.



Western blot testing of human 1) ThP-1, 2) U937, 3) HL60 and 4) HEL cell lysate with AIF-1 antibody. Predicted molecular weight ~17 kDa.