


AIM2 Antibody [clone 812CT9.2.2.1] (F52335)

Catalog No.	Formulation	Size
F52335-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F52335-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgM
Clone Name	812CT9.2.2.1
Purity	Purified
UniProt	O14862
Localization	Cytoplasmic, membrane, nuclear
Applications	Western Blot : 1:1000
Limitations	This AIM2 antibody is available for research use only.



Western blot analysis of lysate from NCI-H292 cell line using AIM2 antibody at 1:1000.
Predicted molecular weight 40-45 kDa.

Description

Tumor suppressor which may act by repressing NF-kappa-B transcriptional activity.

Application Notes

Titration of the AIM2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

This AIM2 antibody was produced from mice immunized with a KLH conjugated synthetic peptide selected from the 10-40 region of human AIM2.

Storage

Aliquot the AIM2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.