


AIM2 Antibody (R34079)


Catalog No.	Formulation	Size
R34079-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9447
Localization	Cytoplasmic, membrane, nuclear
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : 4-6ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This AIM2 antibody is available for research use only.



Western blot testing of IFN β -treated WI-38 lysate with AIM2 antibody at 1ug/ml.
Predicted molecular weight: 40~45kDa.



IHC staining of FFPE human spleen with AIM2 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Absent in melanoma 2

Application Notes

Optimal dilution of the AIM2 antibody should be determined by the researcher.

Immunogen

Amino acids GEKLQLTSGVHSTIK were used as the immunogen for this AIM2 antibody.

Storage

Aliquot and store the AIM2 antibody at -20oC.