


AKAP3 Antibody (R36240)

Catalog No.	Formulation	Size
R36240-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10566
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This AKAP3 antibody is available for research use only.



Western blot testing of human testis lysate with AKAP3 antibody at 1ug/ml. Predicted molecular weight: ~95kDa.

Description

Additional name(s) for this target protein: A-kinase anchor protein 3; sperm oocyte binding protein 1 (SOB1)

Application Notes

Optimal dilution of the AKAP3 antibody should be determined by the researcher.

Immunogen

Amino acids TPLQLLDWLMVNL were used as the immunogen for this AKAP3 antibody.

Storage

Aliquot and store the AKAP3 antibody at -20oC.