


AKAP8 Antibody (R33594)

Catalog No.	Formulation	Size
R33594-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10270
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This AKAP8 antibody is available for research use only.



Western blot testing of human liver lysate with AKAP8 antibody at 1ug/ml. Predicted molecular weight: ~76kDa but routinely observed at ~95kDa.

Description

Additional name(s) for this target protein: A kinase (PRKA) anchor protein 8; A kinase anchor protein 95

Application Notes

Optimal dilution of the AKAP8 antibody should be determined by the researcher.

Immunogen

Amino acids DQGYGGYGAWSAG were used as the immunogen for this AKAP8 antibody.

Storage

Aliquot and store the AKAP8 antibody at -20oC.