


## AKR1C3 Antibody (R33228)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R33228-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                    |                                                               |
|--------------------|---------------------------------------------------------------|
| Availability       | 1-3 business days                                             |
| Species Reactivity | Human                                                         |
| Format             | Antigen affinity purified                                     |
| Clonality          | Polyclonal (goat origin)                                      |
| Isotype            | Goat Ig                                                       |
| Purity             | Antigen affinity                                              |
| Gene ID            | 8644                                                          |
| Applications       | Western Blot : 0.01-0.1ug/ml<br>ELISA (peptide) LOD : 1:32000 |
| Limitations        | This AKR1C3 antibody is available for research use only.      |



Western blot testing of human breast lysate with AKR1C3 antibody at 0.03ug/ml.  
Predicted molecular weight: ~37kDa.

## Description

Additional name(s) for this target protein: Aldo-keto reductase family 1, member C3

## Application Notes

Optimal dilution of the AKR1C3 antibody should be determined by the researcher.

## Immunogen

Amino acids CFASHPNYPYSDEY were used as the immunogen for this AKR1C3 antibody.

## Storage

Aliquot and store the AKR1C3 antibody at -20oC.