


AKT2 Antibody (R35817)

Catalog No.	Formulation	Size
R35817-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	208
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This AKT2 antibody is available for research use only.



Western blot testing of human skeletal muscle lysate with AKT2 antibody at 0.3ug/ml.
Predicted molecular weight: ~56kDa.

Description

Additional name(s) for this target protein: Murine thymoma viral (v-akt) homolog-2

Application Notes

Optimal dilution of the AKT2 antibody should be determined by the researcher.

Immunogen

Amino acids RYDSLGLLELDQRT were used as the immunogen for this AKT2 antibody.

Storage

Aliquot and store the AKT2 antibody at -20oC.