


Alk Antibody [Anaplastic lymphoma kinase] (R36170)

Catalog No.	Formulation	Size
R36170-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	266802
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Alk antibody is available for research use only.



Western blot testing of rat brain lysate with Alk antibody at 0.03ug/ml. Predicted molecular weight ~175kDa.

Description

Additional name(s) for this target protein: Anaplastic lymphoma receptor tyrosine kinase; Anaplastic lymphoma kinase

Application Notes

Optimal dilution of the Alk antibody should be determined by the researcher.

Immunogen

Amino acids TETFHPERLESNS were used as the immunogen for this Alk antibody.

Storage

Aliquot and store the Alk antibody at -20oC.