


ALK1 Antibody (R34483)

Catalog No.	Formulation	Size
R34483-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	94
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This ALK1 antibody is available for research use only.



Western blot testing of MCF7 lysate with ALK1 antibody at 2ug/ml. Predicted molecular weight ~60 kDa.

Description

Additional name(s) for this target protein: Activin receptor-like kinase 1; ACVRL1. This protein is different than the similarly named protein ALK (anaplastic lymphoma kinase), Gene ID 238.

Application Notes

Optimal dilution of the ALK1 antibody should be determined by the researcher.

Immunogen

Amino acids KISNSPEKPKVIQ were used as the immunogen for this ALK1 antibody.

Storage

Aliquot and store the ALK1 antibody at -20oC.