


ALKBH1 Antibody / Nucleic acid dioxygenase ALKBH1 (RQ5956)


Catalog No.	Formulation	Size
RQ5956	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request


Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q13686
Localization	Nuclear, cytoplasmic
Applications	Western Blot : 0.5-1ug/ml Immunofluorescence : 5ug/ml
Limitations	This ALKBH1 antibody is available for research use only.




Immunofluorescent staining of FFPE human U-2 OS cells with ALKBH1 antibody (green) and Alpha Tubulin mAb (red). HIER: steam section in pH6 citrate buffer for 20 min.



Immunofluorescent staining of FFPE human U-2 OS cells with ALKBH1 antibody (green) and DAPI nuclear stain (blue). HIER: steam section in pH6 citrate buffer for 20 min.



Western blot testing of human 1) Jurkat and 2) K562 lysate with ALKBH1 antibody. Predicted molecular weight ~44 kDa.



Western blot testing of human 1) A549, 2) PC-3, 3) Jurkat and 4) K562 cell lysate with ALKBH1 antibody. Predicted molecular weight ~44 kDa.

Description

Nucleic acid dioxygenase ALKBH1 is an enzyme that in humans is encoded by the ALKBH1 gene. It is mapped to 14q24.3. This gene encodes a homolog to the *E. coli* alkB gene product. The *E. coli* alkB protein is part of the adaptive response mechanism of DNA alkylation damage repair. It is involved in damage reversal by oxidative demethylation of 1-methyladenine and 3-methylcytosine.

Application Notes

Optimal dilution of the ALKBH1 antibody should be determined by the researcher.

Immunogen

Amino acids YLK TARVNMTVRQVLATDQNFPLEPIEDEKRD from the human protein were used as the immunogen for the ALKBH1 antibody.

Storage

After reconstitution, the ALKBH1 antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.