


Alpha Actin Antibody (F44155)


Catalog No.	Formulation	Size
F44155-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F44155-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Bovine, Chicken, Mouse, Pig, Rabbit, Rat, Xenopus
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P68133
Applications	Western Blot : 1:1000 IHC (Paraffin) : 1:10-1:50 Immunofluorescence : 1:10-1:50 Flow Cytometry : 1:10-1:50
Limitations	This alpha Actin antibody is available for research use only.




Alpha Actin antibody western blot analysis in (1) A549, (2) RD, and (3) human placenta lysate



Alpha Actin antibody immunohistochemistry analysis in formalin fixed and paraffin embedded human skeletal muscle



Confocal immunofluorescent analysis of alpha Actin antibody with A549 cells followed by Alexa Fluor 488-conjugated goat anti-rabbit IgG (green). DAPI was used as a nuclear counterstain (blue).



Alpha Actin antibody flow cytometric analysis of CEM cells (green) compared to a

