


Alpha Taxilin Antibody / TXLNA (RQ7605)


Catalog No.	Formulation	Size
RQ7605	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request


Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	P40222
Localization	Cytoplasmic, extracellular
Applications	Western Blot : 0.5-1ug/ml Flow Cytometry : 1-3ug/million cells Immunofluorescence : 5ug/ml Direct ELISA : 0.1-0.5ug/ml
Limitations	This Alpha Taxilin antibody is available for research use only.




Immunofluorescent staining of FFPE human PC-3 cells with Alpha Taxilin antibody (red) and DAPI nuclear stain (blue). HIER: steam section in pH6 citrate buffer for 20 min.



Western blot testing of human 1) 293T, 2) K562, 3) HepG2, and 4) HeLa cell lysate with Alpha Taxilin antibody. Predicted molecular weight ~62 kDa but commonly observed at ~75 kDa.



Flow cytometry testing of human U-87 MG cells with Alpha Taxilin antibody at 1ug/million cells (blocked with goat sera); Red=cells alone, Green=isotype control, Blue= Alpha Taxilin antibody.



Flow cytometry testing of human U-251 cells with Alpha Taxilin antibody at 1ug/million cells (blocked with goat sera); Red=cells alone, Green=isotype control, Blue= Alpha Taxilin antibody.

Description

Alpha-taxilin is a protein that in humans is encoded by the TXLNA gene. Predicted to enable syntaxin binding activity. Predicted to be involved in exocytosis. Predicted to act upstream of or within B cell activation. Located in cytoplasm.

Application Notes

Optimal dilution of the Alpha Taxilin antibody should be determined by the researcher.

Immunogen

E. coli-derived recombinant human protein (amino acids M1-A546) was used as the immunogen for the Alpha Taxilin antibody.

Storage

After reconstitution, the Alpha Taxilin antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.