


## Annexin A10 Antibody (R34053)

Catalog No.	Formulation	Size
R34053-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	11199
Applications	Immunohistochemistry (FFPE) : 3-5ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This Annexin A10 antibody is available for research use only.



IHC testing of FFPE human pancreas with Annexin A10 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

### Description

Additional name(s) for this target protein: ANX14, annexin 14; ANXA10

### Application Notes

Optimal dilution of the Annexin A10 antibody should be determined by the researcher.

### Immunogen

Amino acids FCGDYVQGTIFPAPN were used as the immunogen for this Annexin A10 antibody.

## Storage

Aliquot and store the Annexin A10 antibody at -20oC.