


## Annexin A11 Antibody (R36391)


| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R36391-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)


|                    |                                                                                                      |
|--------------------|------------------------------------------------------------------------------------------------------|
| Availability       | 1-3 business days                                                                                    |
| Species Reactivity | Human                                                                                                |
| Format             | Antigen affinity purified                                                                            |
| Clonality          | Polyclonal (goat origin)                                                                             |
| Isotype            | Goat Ig                                                                                              |
| Purity             | Antigen affinity                                                                                     |
| Gene ID            | 311                                                                                                  |
| Applications       | Western Blot : 0.2-1ug/ml<br>Immunohistochemistry (FFPE) : 3-5ug/ml<br>ELISA (peptide) LOD : 1:32000 |
| Limitations        | This Annexin A11 antibody is available for research use only.                                        |



Western blot testing of HeLa lysate with Annexin A11 antibody at 0.2ug/ml. Predicted molecular weight: ~54kDa.



IHC staining of FFPE human placenta with Annexin A11 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.



Western blot testing of A) HeLa and B) Jurkat cell lysate with Annexin A11 antibody.  
Predicted molecular weight: ~54kDa.

## Description

Additional name(s) for this target protein: ANXA11

## Application Notes

Optimal dilution of the Annexin A11 antibody should be determined by the researcher.

## Immunogen

Amino acids YPGYPPPPGGYPP were used as the immunogen for this Annexin A11 antibody.

## Storage

Aliquot and store the Annexin A11 antibody at -20oC.