


Annexin V Antibody (F49707)

Catalog No.	Formulation	Size
F49707-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F49707-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Purified
UniProt	P08758
Applications	Western Blot : 1:1000 IHC (Paraffin) : 1:10-1:50
Limitations	This Annexin V antibody is available for research use only.



Annexin V antibody IHC analysis in formalin fixed and paraffin embedded human hepatocarcinoma.



Western blot analysis of Annexin V antibody and SK-BR-3 lysate. Predicted molecular weight ~36 kDa.

Description

Annexin 5 is an anticoagulant protein that acts as an indirect inhibitor of the thromboplastin-specific complex, which is involved in the blood coagulation cascade. [UniProt]

Application Notes

Titration of the Annexin V antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 262-291 from the human protein was used as the immunogen for this Annexin V antibody.

Storage

Aliquot the Annexin V antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.