


Anti-Histone H3 Antibody (F40721)

Catalog No.	Formulation	Size
F40721-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40721-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Bovine, C. elegans, Chicken, Drosophila, Mouse, Pig, Rat, Xenopus, Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P84243
Applications	Western Blot : 1:1000
Limitations	This anti-Histone H3 antibody is available for research use only.



Description

Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Two molecules of each of the four core histones (H2A, H2B, H3, and H4) form an octamer, around which approximately 146 bp of DNA is wrapped in repeating units, called nucleosomes. The linker histone, H1, interacts with linker DNA between nucleosomes and functions in the compaction of chromatin into higher order structures. The protein targeted by this antibody is a replication-independent member of the histone H3 family.

Application Notes

Titration of the anti-Histone H3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 1-30 from the human protein was used as the immunogen for this anti-Histone H3 antibody.

Storage

Aliquot the anti-Histone H3 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.