


Anti-LAMP2 Antibody (F46338)

Catalog No.	Formulation	Size
F46338-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F46338-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P13473
Applications	Western Blot : 1:1000
Limitations	This anti-LAMP2 antibody is available for research use only.



Western blot analysis of anti-LAMP2 antibody and T47D lysate. The protein can be extensively glycosylated and visualized from 45~110 kDa.

Description

LAMP2 is a member of a family of membrane glycoproteins. This glycoprotein provides selectins with carbohydrate ligands. It may play a role in tumor cell metastasis. It may also function in the protection, maintenance, and adhesion of the lysosome. [Wiki]

Application Notes

Titration of the anti-LAMP2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

This anti-LAMP2 antibody was produced from rabbits immunized with LAMP2 recombinant protein.

Storage

Aliquot the anti-LAMP2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.