


AP2A1 Antibody (Biotin Conjugate) (R33860BTN)

Catalog No.	Formulation	Size
R33860BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat, Pig, Cow
Format	Biotin Conjugate
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q95782
Gene ID	160
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This AP2A1 antibody is available for research use only.



Western blot testing of human frontal cortex lysate with biotinylated AP2A1 antibody at 0.01ug/ml. Predicted molecular weight: ~108 kDa.

Description

Additional name(s) for this target protein: Adaptor-related protein complex 2, alpha 1 subunit; Alpha 1 adaptin

Application Notes

Optimal dilution of the AP2A1 antibody should be determined by the researcher.

1. This AP2A1 antibody is specific for isoform 1 of the protein.

Immunogen

Amino acids ELEPPAPESPMALLADPAPAAD were used as the immunogen for this AP2A1 antibody.

Storage

Aliquot and store the AP2A1 antibody at -20oC.