


APEX2 Antibody (F55013)


Catalog No.	Formulation	Size
F55013-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F55013-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

Bulk quote request


Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity purified
UniProt	Q9UBZ4
Localization	Cytoplasmic, nuclear
Applications	Immunofluorescence : 1:25 Flow Cytometry : 1:10-1:50 (1x10 ⁶ cells) Immunohistochemistry (FFPE) : 1:10-1:50 Western Blot : 1:500-1:1000
Limitations	This APEX2 antibody is available for research use only.




IHC testing of FFPE human lung carcinoma tissue with APEX2 antibody. HIER: steam section in pH6 citrate buffer for 20 min and allow to cool prior to staining.



Immunofluorescent staining of human U-251 cells with APEX2 antibody (green) and anti-Actin (red).



Western blot testing of human MCF7 cell lysate with APEX2 antibody. Predicted molecular weight ~57 kDa.



Flow cytometry testing of human MCF7 cells with APEX2 antibody; Blue=isotype control, Green= APEX2 antibody.

Description

APEX2 may participate in both nuclear and mitochondrial post-replicative base excision repair (BER). In the nucleus functions in the PCNA-dependent BER pathway.

Application Notes

The stated application concentrations are suggested starting points. Titration of the APEX2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 143-171 from the human protein was used as the immunogen for the APEX2 antibody.

Storage

Aliquot the APEX2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.