


APOA1 Antibody (R35834)

Catalog No.	Formulation	Size
R35834-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	335
Localization	Cytoplasmic, membranous
Applications	Immunohistochemistry (FFPE) : 2-4ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This APOA1 antibody is available for research use only.



IHC testing of FFPE human liver with APOA1 antibody at 2ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: Apolipoprotein A-I

Application Notes

Optimal dilution of the APOA1 antibody should be determined by the researcher.

Immunogen

Amino acids SALEEYTKKLNTQ were used as the immunogen for this APOA1 antibody.

Storage

Aliquot and store the APOA1 antibody at -20oC.