


APOA1BP Antibody (R35104)

Catalog No.	Formulation	Size
R35104-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	128240
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This APOA1BP antibody is available for research use only.



Western blot testing of human heart lysate with APOA1BP antibody at 0.5ug/ml.
Predicted molecular weight: ~32kDa.

Description

Additional name(s) for this target protein: Apolipoprotein A-I binding protein

Application Notes

Optimal dilution of the APOA1BP antibody should be determined by the researcher.

Immunogen

Amino acids PPYPDTECVYRLQ were used as the immunogen for this APOA1BP antibody.

Storage

Aliquot and store the APOA1BP antibody at -20oC.