

APOA4 Antibody (R34401)

Catalog No.	Formulation	Size
R34401-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	337
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This APOA4 antibody is available for research use only.



Western blot testing of human duodenum lysate with APOA4 antibody at 0.01ug/ml.
Predicted molecular weight: ~45 kDa.

Description

Additional name(s) for this target protein: Apolipoprotein A-IV

Application Notes

Optimal dilution of the APOA4 antibody should be determined by the researcher.

Immunogen

Amino acids KEKESQDKTSLP were used as the immunogen for this APOA4 antibody.

Storage

Aliquot and store the APOA4 antibody at -20oC.