


APOBEC2 Antibody (R33695)

Catalog No.	Formulation	Size
R33695-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10930
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This APOBEC2 antibody is available for research use only.



Western blot testing of human heart lysate with APOBEC2 antibody at 1ug/ml with [B] and without [A] blocking peptide. Predicted molecular weight: ~26kDa.

Description

Additional name(s) for this target protein: Apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 2

Application Notes

Optimal dilution of the APOBEC2 antibody should be determined by the researcher.

Immunogen

Amino acids EAASQNGEDLENLDD were used as the immunogen for this APOBEC2 antibody.

Storage

Aliquot and store the APOBEC2 antibody at -20oC.