


APOBEC3C Antibody (R35672)

Catalog No.	Formulation	Size
R35672-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	27350
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This APOBEC3C antibody is available for research use only.



Western blot testing of Daudi lysate with APOBEC3C antibody at 0.5ug/ml. Predicted molecular weight: ~23kDa.

Description

Additional name(s) for this target protein: Apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 3C

Application Notes

Optimal dilution of the APOBEC3C antibody should be determined by the researcher.

Immunogen

Amino acids RLLKRRLRESLQ were used as the immunogen for this APOBEC3C antibody.

Storage

Aliquot and store the APOBEC3C antibody at -20oC.