

APOBEC3G Antibody (R33416)


Catalog No.	Formulation	Size
R33416-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	60489
Localization	Cytoplasmic, membranous
Applications	Western Blot : 0.5-1.5ug/ml IHC (FFPE) : 2.5ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This APOBEC3G antibody is available for research use only.



Western blot testing of Daudi lysate with APOBEC3G antibody at 0.5ug/ml. Predicted molecular weight: ~46 kDa.



IHC testing of FFPE human kidney tissue with APOBEC3G antibody at 2.5ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.



IHC testing of FFPE human testis tissue with APOBEC3G antibody at 2.5ug/ml.
Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

Description

Additional name(s) for this target protein: Apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 3G

Application Notes

Optimal dilution of the APOBEC3G antibody should be determined by the researcher.

Immunogen

Amino acids DEHSQDLSGRLR were used as the immunogen for this APOBEC3G antibody.

Storage

Aliquot and store the APOBEC3G antibody at -20oC.