


Apolipoprotein D Antibody (R33804)


Catalog No.	Formulation	Size
R33804-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	347
Applications	Western Blot : 0.3-1ug/ml Immunohistochemistry (FFPE) : 5-10ug/ml Enzyme ImmunoAssay : used as capture for recombinant protein (coat at 1-2ug/ml) ELISA (peptide) LOD : 1:32000
Limitations	This Apolipoprotein D antibody is available for research use only.



Western blot testing of human liver lysate with Apolipoprotein D antibody at 0.3ug/ml. Predicted molecular weight: ~21/21-33kDa (unmodified/glycosylated).



IHC staining of FFPE human kidney with Apolipoprotein D antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: APOD; Apo-D

Application Notes

Optimal dilution of the Apolipoprotein D antibody should be determined by the researcher.

Immunogen

Amino acids EIEKIPTTFENGRC were used as the immunogen for this Apolipoprotein D antibody.

Storage

Aliquot and store the Apolipoprotein D antibody at -20oC.