


Apom Antibody (R34449)

Catalog No.	Formulation	Size
R34449-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	55938
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Apom antibody is available for research use only.



Western blot testing of human placenta lysate with Apom antibody at 0.5ug/ml. Predicted molecular weight: ~21/26kDa (unmodified/glycosylated).

Description

Additional name(s) for this target protein: Apolipoprotein M

Application Notes

Optimal dilution of the Apom antibody should be determined by the researcher.

Immunogen

Amino acids KGNMELRTEGRPD were used as the immunogen for this Apom antibody.

Storage

Aliquot and store the Apom antibody at -20oC.