


APRIL Antibody / ANP32B / PHAPI2 Biotin Conjugate (R34623BTN)

Catalog No.	Formulation	Size
R34623BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat, Pig
Predicted Reactivity	Cow, Dog
Format	Biotin Conjugate
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10541
Applications	Western Blot : 0.1-0.3ug/ml IHC-PARAFFIN PROTOCOL : suitable ELISA (peptide) LOD : 1:16000
Limitations	This APRIL antibody is available for research use only.



Western blot testing of human tonsil lysate with biotinylated APRIL antibody at 0.2ug/ml.
Predicted molecular weight: ~28kDa.

Description

Additional name(s) for this target protein: Acidic protein rich in leucines, SSP29, ANP32B; ANP32B. This protein is NOT the TNF family member also called APRIL (Gene ID 8741).

Application Notes

Optimal dilution of the APRIL antibody should be determined by the researcher.

Immunogen

Amino acids KKKRETDDEGEDD were used as the immunogen for this APRIL antibody.

Storage

Aliquot and store the APRIL antibody at -20oC.