


## ARHGDIG Antibody (R33721)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R33721-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                      |                                                            |
|----------------------|------------------------------------------------------------|
| Availability         | 1-3 business days                                          |
| Species Reactivity   | Human                                                      |
| Predicted Reactivity | Cow, Dog, Mouse, Rat                                       |
| Format               | Antigen affinity purified                                  |
| Clonality            | Polyclonal (goat origin)                                   |
| Isotype              | Goat Ig                                                    |
| Purity               | Antigen affinity                                           |
| Gene ID              | 398                                                        |
| Applications         | Western Blot : 0.5-2ug/ml<br>ELISA (peptide) LOD : 1:16000 |
| Limitations          | This ARHGDIG antibody is available for research use only.  |



Western blot testing of HeLa lysate with ARHGDIG antibody at 0.5ug/ml. Predicted molecular weight: ~25kDa.

## Description

Additional name(s) for this target protein: Rho GDP dissociation inhibitor (GDI) gamma; RHOGDI-3; RhoGDI gamma

## Application Notes

Optimal dilution of the ARHGDIG antibody should be determined by the researcher.

## Immunogen

Amino acids EAVPEYRAPGRK were used as the immunogen for this ARHGDIG antibody.

## Storage

Aliquot and store the ARHGDIG antibody at -20oC.