


ARID5A Antibody (RQ6369)

Catalog No.	Formulation	Size
RQ6369	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q03989
Applications	Western Blot : 0.5-1ug/ml Direct ELISA : 0.1-0.5ug/ml
Limitations	This ARID5A antibody is available for research use only.



Western blot testing of human 1) placenta and 2) K562 cell lysate with ARID5A antibody. Predicted molecular weight ~64 kDa.

Description

Members of the ARID protein family, including ARID5A, have diverse functions but all appear to play important roles in development, tissue-specific gene expression, and regulation of cell growth.

Application Notes

Optimal dilution of the ARID5A antibody should be determined by the researcher.

Immunogen

Recombinant human protein (amino acids R55-L594) was used as the immunogen for the ARID5A antibody.

Storage

After reconstitution, the ARID5A antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.