


## ARK5 Antibody (R35382)

Catalog No.	Formulation	Size
R35382-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9891
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This ARK5 antibody is available for research use only.



Western blot testing of kidney lysate with ARK5 antibody at 0.3ug/ml. Predicted molecular weight: ~75kDa.

## Description

Additional name(s) for this target protein: AMPK-related protein kinase 5; NUAK1

## Application Notes

Optimal dilution of the ARK5 antibody should be determined by the researcher.

## Immunogen

Amino acids QNRPRPQYLKRYRNR were used as the immunogen for this ARK5 antibody.

## Storage

Aliquot and store the ARK5 antibody at -20oC.